

## TRAPEZE BAR PATIENT GUIDELINES

The trapeze bar allows a patient more independence by providing the patient with a means of self help to change position in bed, move onto and off of a bedpan, or transfer into or out of a bed.

Our representative will install the trapeze bar for you, ensuring for tightness of clamps and stability.

*Note:* If any part of the trapeze bar becomes loose or unstable, discontinue use and call our office. Continued use of the trapeze bar may result in serious injury.

Check all clamps daily to ensure that they are securely tight. Adjust the trapeze bar height by hooking the chain on the "S" hook at the correct height. To get the bar out of the way simply loop it over the horizontal bar of the trapeze frame.

### CLEANING INSTRUCTIONS

Clean bar and stand with a damp cloth to keep free from dust.

### FOR MORE INFORMATION

For more information, please call the Carle Medical Supply location nearest you.

1208 N. CUNNINGHAM AVENUE  
URBANA  
(217) 383-3487

1016 BROADWAY  
MATTOON  
(217) 235-6421

622 N. GILBERT  
DANVILLE  
(217) 446-0577



[www.carlemedicalsupply.com](http://www.carlemedicalsupply.com)